

Product Datasheet

Human sCD40L ELISA Kit (orb1658223)

Description	Human sCD40L ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Sensitivity: 15 pg/mL. Gene_name: CD40LG. Gene_full_name: CD40 ligand. Protein_name: CD40 ligand. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: CD40LG. Sequence_similarities: Belongs to the tumor necrosis factor family. Tissue_specificity: Specifically expressed on activated CD4+ T-lymphocytes.
Expiration Date	Please enquire.