

Product Datasheet

Rat PDGF-AB ELISA Kit (orb1658257)

Description	Rat PDGF-AB ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 3 pg/mL. Gene_name: PDGFB PDGFA. Gene_full_name: platelet-derived growth factor beta polypeptide platelet-derived growth factor alpha polypeptide. Protein_name: Platelet-derived growth factor subunit B Platelet-derived growth factor subunit A. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: cardiovascular blood serum proteins platelets signal transduction growth factors hormones pdgf immunology secreted molecules secreted molecules cancer growth factors atherosclerosis vascular inflammation inflammatory mediators neuroscience neurology process neuroregeneration cancer metabolism response to hypoxia pathways and processes metabolism processes types of disease.
Expiration Date	Please enquire.