

Product Datasheet

Rat MIP-3 Alpha/CCL20 ELISA Kit (orb1658275)

Description	Rat MIP-3 Alpha/CCL20 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8 pg/mL-500 pg/mL. Sensitivity: 1 pg/mL. Gene_name: CCL20. Gene_full_name: C-C motif chemokine 20. Protein_name: C-C motif chemokine 20. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. chemokines beta chemokines (cc) signal transduction signaling pathway g protein signaling gpcr kits lysates other kits elisa kits cytokines and cytokine receptors elisa kits. Tissue_specificity: Low levels in thymus and lung.
Expiration Date	Please enquire.