

Product Datasheet

Mouse MIP-2 ELISA Kit (orb1658278)

Description	Mouse MIP-2 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 5 pg/mL. Gene_name: CXCL2. Gene_full_name: C-X-C motif chemokine 2. Protein_name: C-X-C motif chemokine 2. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity chemokines alpha chemokines (cxc).
Expiration Date	Please enquire.