

Product Datasheet

Human CCL4/MIP-1 Beta ELISA Kit (orb1658279)

Description	Human CCL4/MIP-1 Beta ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8 pg/mL-500 pg/mL. Sensitivity: 1 pg/mL. Gene_name: CCL4. Gene_full_name: C-C motif chemokine 4. Protein_name: C-C motif chemokine 4. Sample Type: cell culture supernates serum and plasma (heparin EDTA citrate). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. chemokines beta chemokines (cc) kits lysates other kits elisa kits cytokines and cytokine receptors elisa kits immune system diseases antiviral signaling hiv-related. Sequence_similarities: Belongs to the intercrine beta (chemokine CC) family.
Expiration Date	Please enquire.