

## Product Datasheet

### Human CCL7/MCP-3 ELISA Kit (orb1658284)

<b>Description</b>	Human CCL7/MCP-3 ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 1 pg/mL. Gene_name: CCL7. Gene_full_name: C-C motif chemokine 7. Protein_name: C-C motif chemokine 7. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. chemokines beta chemokines (cc) kits lysates other kits elisa kits cytokines and cytokine receptors elisa kits. Sequence_similarities: Belongs to the intercrine beta (chemokine CC) family.
<b>Expiration Date</b>	Please enquire.