

Product Datasheet

Human IL1R1 ELISA Kit (orb1658304)

Description	Human IL1R1 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 125 pg/mL-8000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: IL1R1. Gene_full_name: Interleukin-1 receptor type 1. Protein_name: Interleukin-1 receptor type 1. Sample Type: cell culture supernates and serum. Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity cytokines interleukins cardiovascular atherosclerosis vascular inflammation inflammatory mediators stem cells hematopoietic progenitors lymphoid t lymphocytic lineage myeloid dendritic cell lineage. Sequence_similarities: Belongs to the interleukin-1 receptor family.
Expiration Date	Please enquire.