

Product Datasheet

Human IGFBP-1 ELISA Kit (orb1658312)

Description	Human IGFBP-1 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 1 pg/mL. Gene_name: IGFBP1. Gene_full_name: Insulin-like growth factor-binding protein 1. Protein_name: Insulin-like growth factor-binding protein 1. Sample Type: cell culture supernates serum plasma (heparin EDTA) and urine. Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: signal transduction growth factors hormones insulin insulin-like cancer growth factors insulin and insulin-like cardiovascular atherosclerosis diabetes associated kits lysates other kits elisa kits growth factors and hormones elisa kits. Sequence_similarities: Contains 1 IGFBP N-terminal domain.
Expiration Date	Please enquire.