

Product Datasheet

Mouse IGF-2 ELISA Kit (orb1658313)

Description	Mouse IGF-2 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.