

Product Datasheet

Human CXCL6/GCP2 ELISA Kit (orb1658321)

Description	Human CXCL6/GCP2 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CXCL6. Gene_full_name: C-X-C motif chemokine 6. Protein_name: C-X-C motif chemokine 6. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: Nothing Found. Sequence_similarities: Belongs to the intercrine alpha (chemokine CxC) family.
Expiration Date	Please enquire.