

Product Datasheet

Human Transmembrane glycoprotein NMB ELISA Kit (orb1658347)

Description	Human Transmembrane glycoprotein NMB ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.15 ng/mL-10.0 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.