

Product Datasheet

CDY1 antibody (orb165872)

Description	Rabbit polyclonal antibody to CDY1 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CDY1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60473
Uniprot ID	Q9Y6F8
NCBI	NP_733841.1 , NP_001003895.1 , NP_004671.1 , NP_001003894.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

