

Product Datasheet

Human Radical S-adenosyl methionine domain-containing protein 2 ELISA Kit (orb1658821)

Description	Human Radical S-adenosyl methionine domain-containing protein 2 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: less than 0.054 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.