

Product Datasheet

Rat Oxytocin receptor ELISA Kit (orb1658993)

Description	Rat Oxytocin receptor ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: less than 0.39 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.