

Product Datasheet

Bovine Prolactin-inducible protein homolog ELISA Kit (orb1659086)

Description	Bovine Prolactin-inducible protein homolog ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 39-2500 pg/mL. Sensitivity: less than 10 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.