

Product Datasheet

Mouse Metastasis-suppressor KiSS-1 ELISA Kit (orb1659167)

Description Mouse Metastasis-suppressor KiSS-1
ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL.
Assay Time: 3-5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
Serum, plasma, tissue homogenates
and other biological fluids.

Expiration Date Please enquire.