

Product Datasheet

NR1H4 antibody (orb165921)

Description	Rabbit polyclonal antibody to NR1H4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NR1H4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55914
Uniprot ID	Q96R11
NCBI	NP_001193906.1 , NP_005114.1 , NP_001193921.1 , NP_001193907.1 , NP_001193908.1 , NP_001193922.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

