

## Product Datasheet

# Human C-C chemokine receptor type 6 ELISA Kit (orb1659216)

<b>Description</b>	Human C-C chemokine receptor type 6 ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.312-20 ng/mL. Sensitivity: less than 0.078 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
<b>Expiration Date</b>	Please enquire.