




Product Datasheet

PMM2 antibody (orb165925)

Description	Rabbit polyclonal antibody to PMM2 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PMM2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	28082
Uniprot ID	O15305
NCBI	NP_000294.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis of HeLa HepG2 K562 MCF-7 cel...



WB analysis of HeLa HepG2 K562 MCF-7 cel...