

Product Datasheet

Human A disintegrin and metalloproteinase with thrombospondin motifs 7 ELISA Kit (orb1659257)

Description	Human A disintegrin and metalloproteinase with thrombospondin motifs 7 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78-5000 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.