

Product Datasheet

Leukadherin-1 (orb16593)

Description	Leukadherin-1 is a specific agonist of CR3 and the leukocyte surface integrin CD11b/CD18.
Note	For research use only
Formula	$C_{22}H_{15}NO_4S_2$
MW	421.4888
CAS Number	[344897-95-6]
SMILES	<chem>S1C(N(C/C1=C([H])/C1=C([H])C([H])=C(C2C([H])=C([H])C(C(=O)O[H])=C([H])C=2[H])O1)=O)C([H])([H])C1C([H])=C([H])C([H])=C([H])C=1[H])=S</chem>
Expiration Date	12 months from date of receipt.