

## Product Datasheet

# Human Tyrosine-protein kinase BTK ELISA Kit (orb1659314)

**Description** Human Tyrosine-protein kinase BTK  
ELISA Kit

**Reactivity** Human

**Tested Applications** ELISA

**Note** For research use only

**Application notes** Detect Range: 0.78-50 ng/mL.  
Sensitivity: less than 0.39 ng/mL.  
Assay Time: 3-5 h. Detection  
wavelength: 450 nm. Sample  
Volume: 100  $\mu$ L. Sample Type:  
Serum, plasma, tissue homogenates  
and other biological fluids.

**Expiration Date** Please enquire.