

Product Datasheet

Human Tubulin beta-3 chain ELISA Kit (orb1659545)

Description	Human Tubulin beta-3 chain ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78-5000 pg/mL. Sensitivity: less than 33 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.