

Product Datasheet

Human Putative testis-specific prion protein ELISA Kit (orb1659795)

Description	Human Putative testis-specific prion protein ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.31-20.0 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.