

Product Datasheet

Human Heat-stable enterotoxin receptor ELISA Kit (orb1659864)

Description	Human Heat-stable enterotoxin receptor ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: less than 0.078 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.