

## Product Datasheet

# Mouse EGF-like module-containing mucin-like hormone receptor-like 1 ELISA Kit (orb1659882)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Mouse EGF-like module-containing mucin-like hormone receptor-like 1 ELISA Kit  |
| <b>Reactivity</b>          | Mouse  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detect Range: 0.156-10 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 $\mu$ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.  |