

## Product Datasheet

# Human Acylamino-acid-releasing enzyme ELISA Kit (orb1660001)

<b>Description</b>	Human Acylamino-acid-releasing enzyme ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.156-10 ng/mL. Sensitivity: less than 0.1 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
<b>Expiration Date</b>	Please enquire.