

Product Datasheet

Mouse IgG receptor FcRn large subunit p51 ELISA Kit (orb1660053)

Description	Mouse IgG receptor FcRn large subunit p51 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.