

## Product Datasheet

# Rat Transient receptor potential cation channel subfamily V member 1 ELISA Kit (orb1660340)

<b>Description</b>	Rat Transient receptor potential cation channel subfamily V member 1 ELISA Kit
<b>Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 78-5000 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 $\mu$ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
<b>Expiration Date</b>	Please enquire.