

Product Datasheet

Rat Peroxisome proliferator-activated receptor gamma coactivator 1-alpha ELISA Kit (orb1660453)

Description	Rat Peroxisome proliferator-activated receptor gamma coactivator 1-alpha ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.