

Product Datasheet

Human Eukaryotic translation initiation factor 2-alpha kinase 3 ELISA Kit (orb1660475)

Description	Human Eukaryotic translation initiation factor 2-alpha kinase 3 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: less than 0.15 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.