

Product Datasheet

Bovine Gamma-interferon-inducible lysosomal thiol reductase ELISA Kit (orb1660503)

| | |
|----------------------------|--|
| Description | Bovine Gamma-interferon-inducible lysosomal thiol reductase ELISA Kit |
| Reactivity | Bovine |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.312-20 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| Expiration Date | Please enquire. |