

Product Datasheet

Human Novel Coronavirus Nucleoprotein (SARS-CoV-2 N) IgG Antibody ELISA Kit (orb1661267)

Description	Human Novel Coronavirus Nucleoprotein (SARS-CoV-2 N) IgG Antibody ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma.
Expiration Date	Please enquire.