

Product Datasheet

Human ABCC11/ATP-binding cassette sub-family C member 11 ELISA Kit (orb1661308)

| | |
|----------------------------|---|
| Description | Human ABCC11/ATP-binding cassette sub-family C member 11 ELISA Kit |
| Reactivity | Human |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 78-5000 pg/mL. Sensitivity: 39.6 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: ABCC11. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| Expiration Date | Please enquire. |