

Product Datasheet

Rat Afm/Afamin ELISA Kit (orb1661317)

Description	Rat Afm/Afamin ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.78-50 ng/mL. Sensitivity: 0.38 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Afm. Sample Volume: 100 μ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.