

## Product Datasheet

# Mouse Adipor1/Adiponectin receptor protein 1 ELISA Kit (orb1661464)

<b>Description</b>	Mouse Adipor1/Adiponectin receptor protein 1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.156 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Adipor1. Sample Volume: 100 $\mu$ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
<b>Expiration Date</b>	Please enquire.