

## Product Datasheet

# Mouse Crhbp/Corticotropin-releasing factor-binding protein ELISA Kit (orb1661482)

<b>Description</b>	Mouse Crhbp/Corticotropin-releasing factor-binding protein ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 6.25-400 ng/mL. Sensitivity: 1.56 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Crhbp. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
<b>Expiration Date</b>	Please enquire.