

Product Datasheet

Rat Adamts7/A disintegrin and metalloproteinase with thrombospondin motifs 7 ELISA Kit (orb1661521)

Description	Rat Adamts7/A disintegrin and metalloproteinase with thrombospondin motifs 7 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.183 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Adamts7. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.