

Product Datasheet

Mouse Clcf1/Cardiotrophin-like cytokine factor 1 ELISA Kit (orb1661629)

Description	Mouse Clcf1/Cardiotrophin-like cytokine factor 1 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 7.8 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Clcf1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.