

## Product Datasheet

# Mouse Ccr1/C-C chemokine receptor type 1 ELISA Kit (orb1661674)

<b>Description</b>	Mouse Ccr1/C-C chemokine receptor type 1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.78-50 ng/mL. Sensitivity: 0.399 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Ccr1. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
<b>Expiration Date</b>	Please enquire.