

## Product Datasheet

# Human ADCYAP1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit (orb1661835)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human ADCYAP1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit  |
| <b>Reactivity</b>          | Human   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detect Range: 6.25-400 pg/mL.<br>Sensitivity: 3.155 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm.<br>Gene_full_name: ADCYAP1.<br>Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.   |