

Product Datasheet

Mouse Agrp/Agouti-related protein ELISA Kit (orb1661851)

| | |
|----------------------------|---|
| Description | Mouse Agrp/Agouti-related protein ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.156-10 ng/mL. Sensitivity: 0.08 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Agrp. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids. |
| Expiration Date | Please enquire. |