

Product Datasheet

Mouse Clu/Clusterin ELISA Kit (orb1661921)

Description	Mouse Clu/Clusterin ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.063 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Clu. Sample Volume: 100 μ L. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.