

Product Datasheet

Human ANXA10/Annexin A10 ELISA Kit (orb1661949)

Description Human ANXA10/Annexin A10 ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL.
Sensitivity: 0.156 ng/mL. Assay Type:
Sandwich. Assay Time: 3-5 h. Detection
wavelength: 450 nm. Gene_full_name:
ANXA10. Sample Volume: 100 µL.
Sample Type: Serum, plasma, tissue
homogenates and other biological
fluids.

Expiration Date Please enquire.