

Product Datasheet

Mouse Cd38/ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 ELISA Kit (orb1662069)

Description	Mouse Cd38/ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 15.68 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Cd38. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.