

Product Datasheet

Rat Cdk5/Cyclin-dependent-like kinase 5 ELISA Kit (orb1662071)

Description	Rat Cdk5/Cyclin-dependent-like kinase 5 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.78-50 ng/mL. Sensitivity: 0.393 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Cdk5. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.