

## Product Datasheet

# Mouse Cebpa/CCAAT/enhancer-binding protein alpha ELISA Kit (orb1662081)

<b>Description</b>	Mouse Cebpa/CCAAT/enhancer-binding protein alpha ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 78-5000 pg/mL. Sensitivity: 39.5 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Cebpa. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
<b>Expiration Date</b>	Please enquire.