

Product Datasheet

Mouse Bid/BH3-interacting domain death agonist ELISA Kit (orb1662107)

Description	Mouse Bid/BH3-interacting domain death agonist ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.78 ng/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Bid. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.