

Product Datasheet

Mouse Ccl20/C-C motif chemokine 20 ELISA Kit (orb1662238)

Description	Mouse Ccl20/C-C motif chemokine 20 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.9-250 pg/mL. Sensitivity: 1.97 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Ccl20. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.