

Product Datasheet

Chicken CCL4/C-C motif chemokine 4 homolog ELISA Kit (orb1662243)

Description Chicken CCL4/C-C motif chemokine
4 homolog ELISA Kit

Reactivity Gallus

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 15.6-1000 pg/mL.
Sensitivity: 5 pg/mL. Assay Type:
Sandwich. Assay Time: 3-5 h.
Detection wavelength: 450 nm.
Gene_full_name: CCL4. Sample
Volume: 100 µL. Sample Type:
Serum, plasma, tissue
homogenates and other biological
fluids.

Expiration Date Please enquire.